


FLOWCHART

MRI CASES FROM HOSPITAL TENGKU AMPUAN AFZAN (HTAA), KUANTAN

1. Before Appointment Day



VERIFIED BY

ASSOC. PROF. DR. AHMAD RAZALI BIN MD RALIB @ MD RAGHIB
MMC No. 92261
Head
Department of Radiology
Sultan Ahmad Shah Medical Centre @IIUM

2. On Appointment Day

Patient will register and get iPesakit registration number at the Registration Counter, SASMEC @IIUM.

Patient will come to Radiology SASMEC @IIUM.

Patient data will be entered and register in RIS

MRI performed

Once completed, the report will be prepared by medical officers and radiologists.

Images will be burn into the CD

Radiology HTAA to come and collect the CD and printed report (every Monday)

Note:

1. No pediatric cases will be accepted, since patient will need sedation or GA.
2. No cases which need sedation or GA will be accepted. If patient has claustrophobia and will need sedation, cases will be canceled.
3. No cases which need comparison with previous MRI or CT will be accepted.
4. No urgent cases will be accepted.
5. Only cold cases will be accepted, unless MRI in HTAA is down and out.

VERIFIED BY



ASSOC. PROF. DR. AHMAD FAZALI BIN MD RALIS @ MD RAGHIEB
MMC No. 32261

Head

Department of Radiology
Sultan Ahmad Shah Medical Centre @IIUM